



POCKET PRINTED FABRICS DESIGN BLACK AND WHITE GEOMETRICS, DE BOISSY

De Boissy
— DEPUIS 1982

Ref. MDPAMOV04



12/10/2024 01:24:31 - For information purposes only, can be subject to change


GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

POCKET PRINTED FABRICS DESIGN BLACK AND WHITE GEOMETRICS, DE BOISSY



Ref. **MDPAMOV04**

Interchangeable pouch and compatible with 3 bag sizes.

Description

It is fixed with a tamper-proof and patented MOLLE fastening system (invented for the US army: MOdular, Lightweight, Loadbearing, Equipment). This ingenious and practical system is well known to firefighters, since it is found on emergency bags.

Made of recycled polyester, GRS certified

-More colorful "jungle" pattern

-All these pouches are interchangeable and compatible with the 3 sizes of bags of the 2 ranges

-They are fixed with a tamper-proof and patented MOLLE fastening system (invented for the US army: MOdular, Lightweight, Loadbearing, Equipment).

This ingenious and practical system is well known to firefighters, since it is found on emergency bags.

-Dimensions w x h x d: 26 x 12 x 1 to 3 cm

- Weight: 0.07 kg

Available in jungle pattern, more colorful.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040219059	HS Code	4202129000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	27	Packed unit / Height (cm)	24
Bare Unit / Height (cm)	12	Packed unit / Length (cm)	28
Bare Unit / Length (cm)	26	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 12/10/2024 01:24:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr