

ISOTHERM POUCH WITH REUSABLE COLD GEL (WITH VELCRO) DE BOISSY

De Boissy
— DEPUIS 1982


Ref. MDPISOT01



 12/10/2024 02:00:37 - For information purposes only, can be subject to change


GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ISOTHERM POUCH WITH REUSABLE COLD GEL (WITH VELCRO) DE BOISSY



Ref. **MDPISOT01**

Description

- Quilted central divider
- Closure with two-way YKK® slide
- 41 elastic loops of different sizes for up to 34 bulbs
- Dimensions w x h x d: 19 x 14 x 4 cm
- Weight: 0.1 kg (with gel bag)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040219783	HS Code	4202129000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.01
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

12/10/2024 02:00:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr