

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_L



14/10/2024 21:17:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_L

Plastic base stool (seat 3.5 cm)

Description

- Seat ø: 43 cm, padded in very high density foam (30 kg / m3)
 - Fireproof seat (fire resistance classification M1)
 - Axis, tube ø: 5.5 cm in stainless steel
 - Black plastic base with 5 branches ø: 58 cm, very stable with double casters
 - Height adjustable by simple pressure and gas piston
 - Height variation from 53 to 70 cm
 - Maximum load: 150kg
- Matching stool to the OSA003 sofa (several colors available for both models)
- Available in Black, Grey and Lavender

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040196138	HS Code	9402900000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Packed unit / Height (cm)	42
Packed unit / Length (cm)	41	Packed unit / Weight	6.2
Bare Unit / Weight (kg)	5.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Lavender	Packed unit / Quantity	1

14/10/2024 21:17:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr