

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_N



14/10/2024 20:41:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_N

Plastic base stool (seat 3.5 cm)

Description

- Seat ø: 43 cm, padded in very high density foam (30 kg / m3)
 - Fireproof seat (fire resistance classification M1)
 - Axis, tube ø: 5.5 cm in stainless steel
 - Black plastic base with 5 branches ø: 58 cm, very stable with double casters
 - Height adjustable by simple pressure and gas piston
 - Height variation from 53 to 70 cm
 - Maximum load: 150kg
- Matching stool to the OSA003 sofa (several colors available for both models)
- Available in Black, Grey and Lavender

Features

| | | | |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040162645 | HS Code | 9402900000 |
| Warranty (months) | 24 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 17 | Packed unit / Height (cm) | 42 |
| Packed unit / Length (cm) | 41 | Packed unit / Weight | 6.2 |
| Bare Unit / Weight (kg) | 5.2 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Black | Packed unit / Quantity | 1 |

14/10/2024 20:41:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr