

## MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507205



12/10/2024 01:52:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507205

=====

### Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446004019	HS Code	90181990
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.398
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Pink	Packed unit / Quantity	1

12/10/2024 01:52:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)