

SPENGLER TB100 GLUCOMETER COMPLETE KIT



Ref. 745000

CE
0537

DISPOSITIF MEDICAL
IVD
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

STERILE



📅 17/10/2024 21:15:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

SPENGLER TB100 GLUCOMETER COMPLETE KIT



Ref. 745000

For self-monitoring and EASY blood glucose management

Description

An easy-to-use device, available in a complete kit.

- Complete measurement in 5 seconds
- Blood sample reduced to 0.5 ul
- 4 programmable alarms (hour + minute)
- Average over 7/14/30 days
- Main site: whole capillary blood at the fingertip
- Alternative sites of possible sampling: on the palm or forearm
- Large LCD screen
- Memories of 300 measurements
- Dimensions: 92 x 58 x 19 mm
- Weight: 60 grams Complete kit consists of:- 1 TB 100

blood glucose meter with CR 2032 3V battery- 10 test strips
 -10 lancets - 1 control solution- 1 lancing pen- 1 kit- 1 user manual, 1

quick user guide, 1 notebook

LPP reimbursement code:
 Glucometer: 6111380
 Strips: 6111350
 Kit: 6111367
 Lancets: 6111373
 Lancing pen: 611139672

Features

CE Marking	CE0537	Medical device	DMDIV (diagnostic In Vitro)
Bar code	3700446041649	HS Code	9027808000
Warranty (months)	36	Sterile	Yes
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	5.8
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	1.9
Packed unit / Length (cm)	16	Bare Unit / Length (cm)	9.2
Packed unit / Weight	0.22	Bare Unit / Weight (kg)	0.6
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	4	Max temperature	30

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr