

Ref. 334129



 20/10/2024 00:34:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)



Ref. 334129

## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

## Features

|                                    |               |                           |            |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking                         | CE            | Medical device            | Class I    |
| Bar code                           | 3700446000011 | HS Code                   | 9018901000 |
| Warranty (months)                  | 6             | Country                   | China      |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |            |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.11       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |            |
| Minimum temperature                | NA            | Max temperature           | NA         |
| Minimum humidity                   | NA            | Maximum humidity          | NA         |
| <b>OTHER FEATURES</b>              |               |                           |            |
| Packed unit / Quantity             | 1             |                           |            |

20/10/2024 00:34:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)