

CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM



Ref. DUCRYO03



 21/10/2024 21:06:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM



Ref. DUCRY003

Effectively treats benign skin lesions

Description

- Cryotherapeutic agent, allowing to deliver a temperature of -55 ° C in less than 15 seconds
- Application time between 15 and 40 seconds
- Complete set including:
 - A 170 ml aerosol can
 - 2mm and / or 5 mm applicator sticks, usable depending on the different lesions
 - An instruction manual
- The sale of this medical equipment is reserved for health institutions, doctors and doctors' offices

Features

CE Marking	CE0344	Medical device	Class IIa
Bar code	8717484005835	HS Code	29091990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	14
Packed unit / Length (cm)	14	Packed unit / Weight	0.3
SPECIFIC STORAGE CONDITIONS			
Max temperature	25		
OTHER FEATURES			

21/10/2024 21:06:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr