


BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-280

CE



21/10/2024 21:06:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-280

=====

Description

- Handle with tip on the screw
- Unscrewable head 2 positions
- Length 25 cm
- Available in blue, black and red

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040215471	HS Code	90189090
Warranty (months)	60	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.089	Bare Unit / Weight (kg)	0.089
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1

21/10/2024 21:06:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr