

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. MEEMER02_N


GARANTIE
2 ANS
2 YEARS
GUARANTEE



 24/10/2024 20:09:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. **MEEMER02_N**

Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment

Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures
- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg
- Available in waterproof version
- Sold empty

Features

| | | | |
|------------------------------------|---------------|-------------------------|--------|
| CE Marking | NONE | Medical device | No |
| Bar code | 8436564542977 | HS Code | 420212 |
| Warranty (months) | 24 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 24 | Bare Unit / Width (cm) | 26 |
| Packed unit / Height (cm) | 22 | Bare Unit / Height (cm) | 19 |
| Packed unit / Length (cm) | 30 | Bare Unit / Length (cm) | 31.6 |
| Packed unit / Weight | 1.3 | Bare Unit / Weight (kg) | 1.29 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Black | Packed unit / Quantity | 1 |

 24/10/2024 20:09:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr