

MINI URITOP READER

Réf. DPTDR50

URITOP+®



CE

DMDIV

MINI URITOP READER

URITOP+®

Réf. DPTDR50

FEATURES

Medical device Yes DMDIV (diagnostic In Vitro)

Marquage CE

Barcode 3661040218045

Guarantee (month): 0

Customs Code 38220010

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) : Width packaging (cm): 21,5

Length (cm) : Length packaging (cm) : 23,5

Height (cm) : Height packaging (cm) : 10,0

Weight (kg) : Weight packaging (kg) : 1,300

SPECIFIC STORAGE CONDITION

Temperature min (N/A) Temperature max (N/A)

Humidity min (N/A) Humidity max (N/A)



21/06/2022 12:30 - This article is available on our website : www.gsh-med.fr with the following reference DPTDR50