

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520630



CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 ANS  
GUARANTEE

📅 24/10/2024 20:29:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520630

=====

### Description

**HIGN ACCURACY**

- Single tubing with connectic with screw
- "Watchmaking" mechanism

**SHOCK RESISTANT**

- Non deformable metal casing
- Shock proof foam manometer protector

**COMFORT**

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

**CUFF**

- Coton or or nylon cuff with velcro fastening
- Sizes: **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm), **M** (26-34 cm), **L** (32-42 cm) ou **XL** (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

### Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007249	HS Code	90189010
Warranty (months)	36	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.122
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

24/10/2024 20:29:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)