

SET OF GUEDEL CANNULAS (NB 1 / 2 / 3)



Ref. 580850



 24/10/2024 20:22:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF GUEDEL CANNULAS (NB 1 / 2 / 3)



Ref. 580850

Description

Set of 3 single-use cannulas in individual sterile bag. cannula for each size n° 1, 2 and 3. The Guedel cannula is used to keep the airway of an unconscious patient open.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446034672	HS Code	902000909009Z
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.035
Bare Unit / Weight (kg)	0.035		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

24/10/2024 20:22:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr