

MINI READER R-50S ANALYSER URITOP

URITOP+®

Ref. DPTDR50



 30/08/2024 15:24:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MINI READER R-50S ANALYSER URITOP

URITOP+®

Ref. DPTDR50

For a quantitative and qualitative reading; and storage of results

Description

- Automatic transfer, reading and interpretation- Compatible with URITOP + 5, 7 and 11 strips- Patient and product traceability- Analysis capacity 40 to 50 samples / hour- Memory 2,000 results- Computer PC connection and waterproof keyboard OPTIONAL: Printing results with the MINI URITOP®

® reader printer sold separately on order

Features

CE Marking	NONE	Medical device	DMDIV (diagnostic In Vitro)
Bar code	3661040218045	HS Code	38220010
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Packed unit / Height (cm)	10
Packed unit / Length (cm)	23	Packed unit / Weight	1.3
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 15:24:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr