

## COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. MEDIA02\_A



 24/10/2024 19:47:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. MEDIA02\_A

Kit designed with the spectrum of patients with diabetes in mind

### Description

- Double insulated compartment to keep the product cool
- Includes a reusable cold gel
- Material: fabric 420D
- Weight: 130 g
- Outdoor measures: 19 x 10 x 5 cm / 0.95 l- Max
- recommended capacity: 600 g
- Available in blue or silver
- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540744	HS Code	4202.12.99.90
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	20	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	14	Bare Unit / Length (cm)	19
Packed unit / Weight	0.25	Bare Unit / Weight (kg)	0.13
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Argente	Packed unit / Quantity	1

24/10/2024 19:47:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)