

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519050



CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
3 ANS  
**3**  
3 YEARS  
GUARANTEE

28/10/2024 20:08:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519050

=====

### Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

### Features

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                           | 3700446024789 | HS Code                   | 90189010                      |
| Warranty (months)                  | 36            | Country                   | France                        |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | 18            | Packed unit / Height (cm) | 8                             |
| Packed unit / Length (cm)          | 26            | Packed unit / Weight      | 1.15                          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | -20           | Max temperature           | 70                            |
| Minimum humidity                   | 85            | Maximum humidity          | 85                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |

28/10/2024 20:08:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)