


NC3 MONITOR NIBP+SPO2 COMEN

Ref. BCMONI01

COMEN



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE



30/10/2024 22:26:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

NC3 MONITOR NIBP+SPO2 COMEN

COMEN

Ref. BCMONI01

The compact, fast and accurate vital signs station, perfect for outpatient at the bedside.

Description

The NC3 monitor measures 2 vital parameters: non-invasive blood pressure (NIBP) and oxygen saturation (SpO2). Particularly fast and easy to use (starts in 1 to 2 seconds, PNI measurement: 20 to 30 seconds) Its weight and small size with optional trolley with wheels or a wall mount make it a very suitable instrument for outpatient at the bedside, in hospitals, clinics or nursing homes, as well as in doctor's private practices.

- 6" LED backlit screen
- PNI and SpO2 accessories for all types of patients (adult, pediatric, neonatal) (delivered with adult PNI cuff and adult SpO2 sensor)
- Measurement of spot check (manual) or programmed (automatic) NIBP
- Integrated rechargeable lithium battery for a continuous autonomy of 12 hours
- Memory of 50 groups / patients data
- Weight: 1.25kg
- Options: Trolley, wall mount


Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219103	HS Code	9018199560
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	13
Packed unit / Length (cm)	31	Bare Unit / Length (cm)	13
Packed unit / Weight	2.95	Bare Unit / Weight (kg)	1.25
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	NA	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

 30/10/2024 22:26:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr