



NC5 MONITOR NIBP+SPO2 COMEN+TEMP

Ref. BCMONI04

COMEN



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3
ANS
3 YEARS
GUARANTEE



30/10/2024 22:46:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

NC5 MONITOR NIBP+SPO2 COMEN+TEMP

COMEN

Ref. BCMONI04

The very intuitive and readable multi-parameter monitor with its 8" touch screen, equipped with a tympanic thermometer

Description

The NC5 monitor is the multi-parameter monitor perfectly suited to all departments in hospitals and clinics, thanks to its dual mode of spot control and bedside monitoring and its small size with optional roller stand or wall mount. - PNI, SpO2 and TEMP measurement (wireless infrared tympanic thermometer)- 8" touch TFT screen for clear display and intuitive use
 -Audible and visual alarms (distinction between patient problem and technical problem, function: nurse call))- PNI and SpO2 accessories for all types of patients (Adult, pediatric and neonatal) (delivered with PNI cuff and SpO2 sensor for adults)
 - Display of early warning scores: NEWS, MEWS, EWS and calculation of clinical parameters
 -Integrated rechargeable 2200mAh Lithium battery with a battery life of 3 hours and autonomy of 8 hours with 4400mAh Lithium battery not integrated
 - Communication: WIFI, HL7, nurse call, with the Hospital's Computer System and the Centralized Monitoring System
 - Options: Foot with wheels, wall mount, shower barcode
 - Weight: 2.3kg

Features

| | | | |
|------------------------------------|---------------|-------------------------|------------|
| CE Marking | CE1639 | Medical device | Class IIb |
| Bar code | 3661040219080 | HS Code | 9018199560 |
| Warranty (months) | 36 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 29 | Bare Unit / Width (cm) | 16.5 |
| Packed unit / Height (cm) | 26 | Bare Unit / Height (cm) | 25 |
| Packed unit / Length (cm) | 36 | Bare Unit / Length (cm) | 16.5 |
| Packed unit / Weight | 3.8 | Bare Unit / Weight (kg) | 2.3 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 55 |
| Minimum humidity | NA | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

 30/10/2024 22:46:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr