


NC5 MONITOR NIBP+SPO2 COMEN+TEMP+ECG+PRINT

Ref. BCMONI05

COMEN



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3
ANS
3 YEARS
GUARANTEE



30/10/2024 22:31:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

NC5 MONITOR NIBP+SPO2 COMEN+TEMP+ECG+PRINT

COMEN

Ref. BCMONI05

The multi-parameter monitor is very intuitive and highly readable thanks to its 8" touch screen, equipped with a tympanic thermometer and an ECG with printer.

Description

The NC5 monitor is the multi-parameter monitor perfectly suited to all departments in hospitals and clinics, thanks to its dual mode of spot control and bedside monitoring and its small size with optional roller stand or wall mount.

- PNI, SpO2, TEMP, ECG + Printer measurement
- 8" touch TFT screen for clear display and intuitive operation
- Body temperature measurement by wireless infrared tympanic thermometer
- ECG measurement 3 leads with printer
- Audible and visual alarms (distinction between patient problem and technical problem, nurse call function)
- PNI and SpO2 accessories for all types of patients (Adult, pediatric and neonatal) (comes with PNI armband and SpO2 sensor for adults)
- Integrated rechargeable 2200mAh Lithium battery with a battery life of 3 hours and an autonomy of 8 hours with 4400mAh Lithium battery not integrated.
- Communication: WIFI, HL7, nurse call with the Hospital Computer System and the Centralized Monitoring System
- Options: Foot with wheels, wall mount and shower barcode
- Weight: 2,3kg

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219097	HS Code	9018199560
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	29	Bare Unit / Width (cm)	16.5
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	36	Bare Unit / Length (cm)	16.5
Packed unit / Weight	3.8	Bare Unit / Weight (kg)	2.3
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 30/10/2024 22:31:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr