



EMERGENCY/TRANSPORT MONITOR C30 STANDARD CONFIGURATION COMEN

Ref. BCMONI07



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3
ANS
3 YEARS
GUARANTEE



30/10/2024 22:49:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

EMERGENCY/TRANSPORT MONITOR C30 STANDARD CONFIGURATION COMEN

COMEN

Ref. BCMONI07

A 5-in-1 monitor with monitoring of 5 vital parameters. Intuitive thanks to its color touch screen and very easy to carry because ultra compact.

Description

The C30 multi-parameter monitor, dedicated to first aid for on-site monitoring or in emergency vehicles. It measures non-invasive blood pressure (PNI), oxygen saturation (SpO2), respiratory frequency (RESP), capnography (EtCO2), body temperature and has a 3-lead ECG.

Parameters: PNI, SpO2, ECG, Respiration, 2 temperature channels (probes to be added), EtCO2 (sensor to be added).

Screen: 4.3" touch with backlight.

Double alarm: audible (adjustable) and 3-level visual with 3-color smart alarm system (blue, yellow, red).

PNI: suitable for all types of patients (adult, child and newborn).

Storage: 48 hours of data storage and trend analysis.

Battery: Built-in 2200mAh lithium (4 hours battery life).

Options: Optional carrying bag and carrying battery kit (7200mAh). Power

supply: Power kit (12V) to emergency vehicle, optional.

Operating system: Linux.

Dimensions: 19 x 8.2 x 10.5cm.

Weight: 1.5kg.

Warranty: 3 years.

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219479	HS Code	9018193010
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	26	Bare Unit / Width (cm)	8.2
Packed unit / Height (cm)	21	Bare Unit / Height (cm)	10.5
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	19
Packed unit / Weight	3.2	Bare Unit / Weight (kg)	1.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	60
Minimum humidity	NA	Maximum humidity	93
OTHER FEATURES			
Packed unit / Quantity	1		

 30/10/2024 22:49:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr