


NC5 MONITOR NIBP + SPO2 COMEN + TEMP + PRINT

COMEN

Ref. BCMONI09



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3
ANS
3 YEARS
GUARANTEE



30/10/2024 22:40:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

NC5 MONITOR NIBP + SPO2 COMEN + TEMP + PRINT

COMEN

Ref. BCMONI09

The very intuitive and readable multi-parameter monitor with its 8" touch screen, equipped with an eardrum thermometer and a printer

Description

The NC5 monitor is the multi-parameter monitor perfectly suited to all departments of hospitals and clinics, thanks to its dual mode of spot control and continuous monitoring at the patient's bedside and its small size with optional roller foot or wall mount.

- PNI, SpO2 and TEMP (wireless infrared tympanic thermometer) measurement.
- prints the results.
- 8" touch TFT screen for clear display and intuitive use
- Audible and visual alarms (distinction between physiological and technical problems, function: nurse call).
- PNI and SpO2 accessories for all types of patients (Adult, pediatric and neonatal) (comes with PNI cuff and SpO2 sensor for adults).
- Display of early warning scores: NEWS, MEWS, EWS and calculation of clinical parameters.
- Built-in rechargeable 2200mAh lithium battery with 3 hours of battery life and 8 hours of battery with 4400mAh lithium battery not integrated.
- Communication: WIFI, HL7, nurse call with Hospital Computer System (HIS) and Centralized Monitoring System.
- Options: Foot with wheels, wall mount, shower barcode
- Weight: 2.3kg

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040220017	HS Code	9018199560
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	29	Bare Unit / Width (cm)	16.5
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	36	Bare Unit / Length (cm)	16.5
Packed unit / Weight	3.8	Bare Unit / Weight (kg)	2.3
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 30/10/2024 22:40:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr