

## FEMALE TUNIC, IRIS LYOCELL , S0, HOLTEX

Ref. WMIRI07\_96



30/10/2024 22:33:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## FEMALE TUNIC, IRIS LYOCELL , S0, HOLTEX



Ref. WMIRI07\_96

### The eco-responsible women's medical tunic

### Description

- White women's tunic with short sleeves,
- Eco-responsible: 50% Lyocell (eucalyptus fiber) and 50% Polyester,
- Rounded neckline,
- Side opening with snap buttoning
- Nice finishes with white edging with 3 plated pockets and Kimono sleeves,
- Twill, 215g/m2,
- Available from Size 0 to Size 6.

The Holtex eco-friendly and responsible medical tunic brings extra comfort and strength to a classic poly-cotton tunic: superior quality, increased strength, easy maintenance and maximum ease. Its feminine cut is traditional and modern at the same time, perfectly corresponds to the expectations of medical and paramedical professionals: nurses, nursing assistants, doctors, nurses, physiotherapists, dental assistants, audioprosthodontists, dieticians, podiatrists, speech therapists, osteopaths, etc.

Lyocell is a natural material, made from eucalyptus wood pulp.

- Soft, lightweight, comfortable and slightly elastic,
- Easy to maintain, wrinkle resistant, it does not lint.
- Antibacterial, hypoallergenic, breathable, it absorbs moisture (up to 50% more than cotton): it is ideal for sensitive skin.
- Environmentally friendly: Eucalyptus cultivation requires much less water and pesticides than cotton.

### Features

CE Marking	NONE	Medical device	No
Bar code	3661041009468	HS Code	6211433100
Country	Tunisia		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Bleaching	Bleaching prohibited	Washing	Wash at 60°C
Model	Female model	Cleaning	Dry cleaning
Packed unit / Quantity	1	Size	Size 6

30/10/2024 22:33:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)