

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.112



CE



04/11/2024 23:12:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED[®]

Ref. M1.111.112

Reflex Hammer according to Tromner | large

Description

- Reflex hammer made of aluminium / brass
-

Optimal weight distribution

- With neurological sensor tip
- Length: 23 cm | weight: 245 g
- Without brush due to hygienic reasons

Features

| | | | |
|------------------------------------|----------------|---------------------------|---------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381960287 | Warranty (months) | 24 |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -10 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Colors | Black | | |

04/11/2024 23:12:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr