

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.915



04/11/2024 20:37:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX HAMMER ACCORDING TO TROMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.915

Reflex Hammer according to Tromner | large

Description

- Reflex hammer made of aluminium / brass
-

Optimal weight distribution

- With neurological sensor tip
- Length: 23 cm | weight: 245 g
- Without brush due to hygienic reasons

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960317	Warranty (months)	24
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	White		

📅 04/11/2024 20:37:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr