

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330002



 04/11/2024 22:34:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330002

Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446000264	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.274
Bare Unit / Weight (kg)	0.274		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

04/11/2024 22:34:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr