

## PULSE STETHOSCOPE / RED

Ref. 507500



CE



04/11/2024 21:59:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)


+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PULSE STETHOSCOPE / RED



Ref. 507500



### Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009267	HS Code	90181990
Warranty (months)	36		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	24	Packed unit / Weight	0.22
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

04/11/2024 21:59:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)