

TWILL CHOCKER, WHITE, SHORT SLEEVED, SIZE 00


Ref. WSMED07_00



 04/11/2024 20:55:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TWILL CHOCKER, WHITE, SHORT SLEEVED, SIZE 00



Ref. **WSMED07_00**

Description

- Women's tunic- Twill
- Round neckline
- Short kimono sleeves
- 3 plated pockets
- Offset closure on the front by buttoning, 5 pressures
- Height clothing +/- 75 cm
- 65% polyester, 35% cotton
- 200 g / m²
- Available from size 00 to size 6

Features

CE Marking	NONE	Medical device	No
Bar code	3661040197883	HS Code	62114310030
Country	Tunisia		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Packed unit / Height (cm)	1
Packed unit / Length (cm)	38	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.3		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Bleaching	Bleaching prohibited	Washing	Wash at 80°C
Model	Female model	Cleaning	Dry cleaning
Packed unit / Quantity	1	Size	Size 00

📅 04/11/2024 20:55:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr