

## CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)



Ref. 260051



13/11/2024 00:16:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)



Ref. 260051

Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 10" tablet, a product of a new technology, totally autonomous.

### Description

A very great simplicity of use: - 10 "color touch screen- Portable: compact size and light weight, easy to carry: 1.2 kg and 272.6 x 181.2 x 25.8 mm)- Large capacity battery, rechargeable on mains or on the recording base- Allows to print ECGs directly on a printer via wifi network or on the recorder base High-performance technology:

High accuracy: sampling speed 32,000 Hz

- Detection and marking of the pacemaker

- Automatic interpretation (Glasgow) and calculation of intervals- Input of patient information using a keyboard or barcode reader-

Visualization of the 12 tracks (15/18 other tablet optional) at the same time as automatic interpretation and interval calculations.

- Data export via wifi or USB port under .pdf or .xml

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6971049341796	HS Code	9018110000
Warranty (months)	60	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	24	Packed unit / Height (cm)	36
Packed unit / Length (cm)	36	Packed unit / Weight	7
Bare Unit / Weight (kg)	3.3		

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

#### OTHER FEATURES

📅 13/11/2024 00:16:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

🌐 www.gsh-med.fr