

Ref. C1.716.114



Technologie CCT  
2 couleurs de lumière

+ Graduation



12/11/2024 23:20:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

Ref. C1.716.114

### Description

### Features

CE Marking	NONE	Medical device	No
Bar code	04260381963622	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	Black		

12/11/2024 23:20:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)