

## CARDIOMATE 6 - ELECTROCARDIOGRAPH 6 CHANNELS



Ref. 260006



13/11/2024 00:12:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# CARDIOMATE 6 - ELECTROCARDIOGRAPH 6 CHANNELS



Ref. 260006

**CARDIOMATE® ECGs are highly reliable instruments.**

## Description

ECG Cardiomate 3® and 6® are robust and reliable ECG.

Featuring a 7" color screen tilt they allow viewing 12 lead simultaneously.

Simple to use, recording and printing ECG are triggered with only one button.

The Cardiomate 3® & 6® are Cardiomate equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer.

The Cardiomate 3® & Cardiomate 6® are equipped with automatic interpretation and pacemaker detection.

The internal battery life of 3 hours, their lightness and carrying handle make Cardiomate 3® & Cardiomate 6® the perfect companion in the office or traveling.

- Carrying handheld.
- Tilting screen : 7-inch color LCD display (resolution: 800x480 px)
- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM
- Autonomy :The operating life is 1080 ECG.
- Automatic interpretation and intervals calculation.
- Pacemaker detection.
- Network connection, USB & printer.
- Printing and rapid analysis.
- Acquisitions and display 12 simultaneous tracks.
- Alphanumeric keypad.

## Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446025489	HS Code	9018110000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	24.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	31
Packed unit / Weight	4.86	Bare Unit / Weight (kg)	2.4
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	50
Minimum humidity	25	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 13/11/2024 00:12:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)