

## ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99306



12/11/2024 23:08:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. **DS99306**

### Description

Adult cuff (26 to 35 cm)

Washable nylon cuff without pocket with Velcro fastening

### Features

CE Marking	CE0483	Medical device	Class Im (measuring function)
Bar code	3661040075044	HS Code	9018901000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.093	Bare Unit / Weight (kg)	0.093
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

12/11/2024 23:08:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)