

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®

Ref. M1.121.112



12/11/2024 23:28:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED[®]

Ref. M1.121.112

Reflex Hammer according to Buck | small

Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960249	Warranty (months)	24
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		

 12/11/2024 23:28:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr