

## GP DOCTOR BAG BLACK



Ref. MEGP03\_B



12/11/2024 22:49:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## GP DOCTOR BAG BLACK



Ref. **MEGP03\_B**

### Description

Lightweight and easy to use bag thanks to its internal velcro dividers that allow to organize the sanitary equipment in an easy and fast way. Its system of integrated shoulder straps makes it possible to transform it into a backpack.

Light and resistant textile material, anti-tear (Ripstop)

- Special area for personal belongings, business cards, pen, etc ...

- Document compartment, A4

- Internal velcro dividers

- Insulated amp maker, 43 bulbs

- Zipper with eyelet for padlock.

- Large external pocket with double compartment and elastic bands

- Integrated shoulder straps for transformation into a backpack

- Weight: 1 kg

- Outdoor measures: 40 x 21 x 25 cm / 21 L

- Maximum recommended capacity: 3.5 Kg

- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540454	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	40
Packed unit / Weight	1.5	Bare Unit / Weight (kg)	1

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Colors	Gris/turquoise
--------	----------------

12/11/2024 22:49:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

[www.gsh-med.fr](http://www.gsh-med.fr)