

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. DS99606



13/11/2024 00:31:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. **DS99606**

Description

Adult cuff (26 to 35 cm)

Washable nylon cuff without pocket with Velcro fastening

Features

CE Marking	CE0483	Medical device	Class Im (measuring function)
Bar code	3661040075310	HS Code	9018901000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.8	Bare Unit / Weight (kg)	0.122
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/11/2024 00:31:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr