

Ref. 334111




CE



 14/11/2024 22:02:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr



Ref. 334111

Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446025052	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

14/11/2024 22:02:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr