

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519001



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

14/11/2024 22:18:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519001

=====

Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007225	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.905
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 14/11/2024 22:18:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr