

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519029



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

14/11/2024 21:47:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519029

=====

Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

Features

| | | | |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking | CE0459 | Medical device | Class Im (measuring function) |
| Bar code | 3700446008734 | HS Code | 90189010 |
| Warranty (months) | 36 | Country | France |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.663 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 70 |
| Minimum humidity | 85 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

14/11/2024 21:47:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr