

## OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

### LUXAMED®

Ref. A1.426.114




Rechargeable USB



14/11/2024 21:33:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®

Ref. A1.426.114

Innovative otoscope in pocket size

### Description

Illumination: LED-Ring with 6 LEDs (concentric),

no fibre optic, field of vision is not restricted

- Maintenance-free, because no lamp replacement necessary
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours

by constant current drive

- Intensity of illumination: approx. 10.000 Lux
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Pivoting 3-fold magnifying glass with locking function
- Insufflation port for pneumatic test
- Two-component handle (aluminium / fibre-glass reinforced),

also combinable with LuxaScope LED dermatoscope- and

ophthalmoscope head

- Electronic ON/OFF push-button switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm
- Incl. recharger with USB plug (adapter EU/UK/US)

 14/11/2024 21:33:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

**LUXAMED®**

Ref. A1.426.114


### Features

|                                    |                |                           |          |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking                         | CE             | Medical device            | Class I  |
| Bar code                           | 04260381960041 | HS Code                   | 90189084 |
| Warranty (months)                  | 24             |                           |          |
| <b>DIMENSIONS AND WEIGHT</b>       |                |                           |          |
| Packed unit / Width (cm)           | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc             | Packed unit / Weight      | nc       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |                |                           |          |
| Minimum temperature                | 5              | Max temperature           | 40       |
| Minimum humidity                   | 10             | Maximum humidity          | 90       |
| <b>OTHER FEATURES</b>              |                |                           |          |
| Colors                             | Black          |                           |          |

 14/11/2024 21:33:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)