

OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®

Ref. A1.426.214



Rechargeable USB



14/11/2024 22:25:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®

Ref. A1.426.214

Innovative otoscope in pocket size

Description

Illumination: LED-Ring with 6 LEDs (concentric),

no fibre optic, field of vision is not restricted

- Maintenance-free, because no lamp replacement necessary
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours

by constant current drive

- Intensity of illumination: approx. 10.000 Lux
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Pivoting 3-fold magnifying glass with locking function
- Insufflation port for pneumatic test
- Two-component handle (aluminium / fibre-glass reinforced),

also combinable with LuxaScope LED dermatoscope- and

ophthalmoscope head

- Electronic ON/OFF push-button switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm
- Incl. recharger with USB plug (adapter EU/UK/US)

 14/11/2024 22:25:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®

Ref. **A1.426.214**

Features

| | | | |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381960058 | HS Code | 90189084 |
| Warranty (months) | 24 | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | 5 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Colors | Blue | | |

 14/11/2024 22:25:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr