

Ref. C1.616.914



Technologie CCT  
2 couleurs de lumière



14/11/2024 21:30:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

Ref. C1.616.914

Innovative dermatoscope in pocket size

### Description

Innovative illumination: LED-Ring with 2x5 LEDs in concentric

structure for a highly detailed image, homogenous light, switchable

with one button press : warm white light or cold white light


- Removable contact plate for easy cleaning
- Maintenance-free, because no lamp replacement necessary
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
- Two-component-handle (aluminium/fibre-glass reinforced),
- Also combinable with LuxaScope Auris CCT LED 2.5 V

otoscope-head

14/11/2024 21:30:02 - For information purposes only, can be subject to change  
• Electronic ON/OFF push-button-switch

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

Ref. C1.616.914

### Features

CE Marking	NONE	Medical device	No
Bar code	04260381963837	HS Code	90189084
Warranty (months)	24		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	White		

14/11/2024 21:30:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)