

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED**<sup>®</sup>

Ref. D1.212.412

CE



 14/11/2024 21:17:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**

Ref. **D1.212.412**

Penlight for eye and throat diagnostic., long life technology.

### Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

### Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961772	Warranty (months)	24
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	Red		

14/11/2024 21:17:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)