

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. D1.312.912

CE



 14/11/2024 21:46:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. **D1.312.912**

Penlight for eye and throat diagnostic., long life technology.

Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Features

| | | | |
|------------------------------------|----------------|---------------------------|---------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381963349 | Warranty (months) | 24 |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -10 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Colors | Kids | | |

14/11/2024 21:46:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr