

CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E



Ref. DPOXY04



 14/11/2024 21:28:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E



Ref. DPOXY04

Description

With its remote pediatric sensor, flexible, reusable and immersible, it will adapt to children from 2 to 8 years old. The PC 60E pulse oximeter is reliable and accurate. It for health professionals in the city.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040163833	HS Code	9018909911
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	2
Packed unit / Length (cm)	14	Packed unit / Weight	0.001
Bare Unit / Weight (kg)	0.059		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 14/11/2024 21:28:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr