

BLACK RING FOR MAGISTER/PULSE BELL



Ref. 507941



 26/11/2024 09:40:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BLACK RING FOR MAGISTER/PULSE BELL



Ref. 507941

Description

Spare part for Pulse® & Dual Pulse® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009700	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/11/2024 09:40:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr