

## PHIAL AMPOULARIUM ELITE BAGS EB09.002



Ref. MEAMPO01



27/08/2024 08:30:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PHIAL AMPOULARIUM ELITE BAGS EB09.002



Ref. MEAMPO01

Insulated amp maker, 43 bulbs

### Description

- Padded central separator  
-

Fasteners- Color: black

- Weight: 160 g  
-

Outdoor measures: 18 x 13 x 4 cm / 0.93 L  
-

Maximum recommended capacity: 1 kg

- Sold empty

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540546	HS Code	420212
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	13
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	18	Bare Unit / Length (cm)	18
Packed unit / Weight	0.276	Bare Unit / Weight (kg)	0.16
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 08:30:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)