

FINSTERER SUCTION TUBE 32 CM

Ref. IC46129



26/08/2024 22:32:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

FINSTERER SUCTION TUBE 32 CM



Ref. IC46129

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040011790	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	32.5
Packed unit / Weight	0.136	Bare Unit / Weight (kg)	0.132
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 22:32:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr