

ASKIR C30 ASKIR C30 SUCTION UNIT 2X2L



Ref. BCASPI33



29/08/2024 14:50:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ASKIR C30 ASKIR C30 SUCTION UNIT 2X2L



Ref. BCASPI33

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Description

Ca-mi Askir C30 Suction Unit, autoclavable polycarbonate jars, 2 x 2000 ml

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Intended for professional use and equipped with a powerful suction system that can drain up to 40 l/min max, Askir C30 is equipped with two collection jars and a five-caster system making it easy to transport.

Available in Basic or FS

version - BASIC version, including: Fluid collection jar with anti-reflux valve / 8 x 14 Silicone hoses (autoclavable) and tapered connectors ø 10-11-12 mm / Antibacterial & hydrophobic filter (single patient) / Power cord with Schuko plug -

Main system equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 40 l/min

-

Noise level: 61 dBA

- Weight: 6.2 kg / Dimensions: 32 x 99 x 30 cm

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139197	HS Code	90189084

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

+33(0)4 42 90 31 31

info@gsh-med.fr

www.gsh-med.fr

Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			