

LITHIUM BATTERY

Ref. SJACCE32

SaverOne®


CE



 29/08/2024 13:15:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LITHIUM BATTERY



Ref. **SJACCE32**

Description

For Défibrilateur Saver One, new generation, ie after November 2013.

- Capacity: 3500 Ah
- Battery voltage: 25.2 V
- Non-rechargeable battery

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040115450	HS Code	85076000
Warranty (months)	48	Country	Hungary
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Bare Unit / Width (cm)	7.7
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3.1
Packed unit / Length (cm)	8	Bare Unit / Length (cm)	21.7
Packed unit / Weight	0.35	Bare Unit / Weight (kg)	0.3
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 13:15:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr