

KERRISSON PUNCH FORCEPS UP 90° 20 CM X 4 MM

Ref. IP46841



27/08/2024 04:36:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

KERRISSON PUNCH FORCEPS UP 90° 20 CM X 4 MM



Ref. IP46841

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040026411	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	9.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.182	Bare Unit / Weight (kg)	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 04:36:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr