

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 180 MM - 2 MM



Ref. IP46846



27/08/2024 04:38:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 180 MM - 2 MM



Ref. IP46846

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040026466	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	1
Packed unit / Length (cm)	27	Packed unit / Weight	0.208
Bare Unit / Weight (kg)	0.208		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 04:38:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr